

FXR1 Antibody / Fragile X Mental Retardation 1 (R35634)

Catalog No.	Formulation	Size
R35634-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Rat, Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	8087
Applications	Western Blot : 0.2-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This FXR1 antibody is available for research use only.



Western blot testing of NIH3T3 cell lysate with FXR1 antibody at 0.2ug/ml.
Predicted/observed molecular weight ~60/80kDa.

Description

Additional name(s) for this target protein: Fragile X mental retardation-related protein 1; FMRP; FMR1

Application Notes

Optimal dilution of the FXR1 antibody should be determined by the researcher.

Immunogen

Amino acids RIEGDNENKLPRED were used as the immunogen for this FXR1 antibody.

Storage

Aliquot and store the FXR1 antibody at -20oC.