


FYB Antibody (R35871)

Catalog No.	Formulation	Size
R35871-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2533
Applications	Western Blot : 0.2-0.6ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This FYB antibody is available for research use only.



Western blot testing of MOLT4 lysate with FYB antibody at 0.2ug/ml. Predicted molecular weight: ~90kDa.

Description

Additional name(s) for this target protein: FYN-binding protein

Application Notes

Optimal dilution of the FYB antibody should be determined by the researcher.

Immunogen

Amino acids SEEKPDKEPKP were used as the immunogen for this FYB antibody.

Storage

Aliquot and store the FYB antibody at -20oC.