


FYN Antibody / Tyrosine-protein kinase Fyn (R36336)

Catalog No.	Formulation	Size
R36336-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Rat, Dog, Pig, Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2534
Applications	Western Blot : 0.2-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This FYN antibody is available for research use only.



Western blot testing of A431 [A] and Jurkat lysate [B] with FYN antibody at 0.3ug/ml.
Predicted molecular weight: ~59/59/54 kDa (isoforms 1/2/3).

Description

Additional name(s) for this target protein: FYN oncogene related to SRC FGR YES; src-like kinase (SLK)

Application Notes

Optimal dilution of the FYN antibody should be determined by the researcher.

Immunogen

Amino acids VQCKDKEATKLTE were used as the immunogen for this FYN antibody.

Storage

Aliquot and store the FYN antibody at -20oC.