


## FYN Antibody / Tyrosine-protein kinase Fyn (R36336)

Catalog No.	Formulation	Size
R36336-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human, Mouse
<b>Predicted Reactivity</b>	Rat, Dog, Pig, Cow
<b>Format</b>	Antigen affinity purified
<b>Host</b>	Goat
<b>Clonality</b>	Polyclonal (goat origin)
<b>Isotype</b>	Goat Ig
<b>Purity</b>	Antigen affinity
<b>Gene ID</b>	2534
<b>Applications</b>	Western Blot : 0.2-1ug/ml ELISA (peptide) LOD : 1:32000
<b>Limitations</b>	This FYN antibody is available for research use only.



Western blot testing of A431 [A] and Jurkat lysate [B] with FYN antibody at 0.3ug/ml.  
Predicted molecular weight: ~59/59/54 kDa (isoforms 1/2/3).

### Description

Additional name(s) for this target protein: FYN oncogene related to SRC FGR YES; src-like kinase (SLK)

### Application Notes

Optimal dilution of the FYN antibody should be determined by the researcher.

## **Immunogen**

Amino acids VQCKDKEATKLTE were used as the immunogen for this FYN antibody.

## **Storage**

Aliquot and store the FYN antibody at -20oC.