


## FZD4 Antibody (R34482)

Catalog No.	Formulation	Size
R34482-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	8322
Applications	Immunohistochemistry (FFPE) : 2-4ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This FZD4 antibody is available for research use only.



IHC testing of FFPE human kidney with FZD4 antibody at 2ug/ml. HIER: steamed with pH6 citrate buffer, HRP-staining.

## Description

Additional name(s) for this target protein: Frizzled 4

## Application Notes

Optimal dilution of the FZD4 antibody should be determined by the researcher.

## Immunogen

Amino acids KIRSNLQKDGTKT were used as the immunogen for this FZD4 antibody.

## Storage

Aliquot and store the FZD4 antibody at -20oC.