


GAB2 Antibody (R35496)

Catalog No.	Formulation	Size
R35496-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9846
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This GAB2 antibody is available for research use only.



Western blot testing of K562 cell lysate with GAB2 antibody at 0.5ug/ml. Predicted molecular weight: ~75kDa.

Description

Additional name(s) for this target protein: GRB2-associated binding protein 2

Application Notes

Optimal dilution of the GAB2 antibody should be determined by the researcher.

Immunogen

Amino acids QVDKEKTQALQN were used as the immunogen for this GAB2 antibody.

Storage

Aliquot and store the GAB2 antibody at -20oC.