

Zebrafish Gabarapa Antibody (F53276)

Catalog No.	Formulation	Size
F53276-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F53276-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Zebrafish
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q6PSS4
Applications	Western Blot : 1:500
Limitations	This Gabarapa antibody is available for research use only.



Western blot testing of Gabarapa antibody at 1:500 dilution + zebrafish lysate; Predicted molecular weight : 14 kDa.

Description

Zebrafish Gabarapa Antibody detects GABA(A) receptor associated protein.

This Zebrafish antibody is part of a [broader Zebrafish / Danio rerio antibody panel](#) offered by NSJ Bioreagents.

Application Notes

Titration of the Gabarapa antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

This Gabarapa antibody was produced from a rabbit immunized with a KLH conjugated synthetic peptide between 33-66 amino acids from the human region of zebrafish Gabarapa.

Storage

Aliquot the Gabarapa antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.