


GABPB2 Antibody (R36065)


Catalog No.	Formulation	Size
R36065-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)


Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2553
Localization	Cytoplasmic, nuclear
Applications	Western Blot : 0.5-1.5ug/ml IHC (FFPE) : 3.75ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This GABPB2 antibody is available for research use only.



Western blot testing of HeLa cell lysate with GABPB2 antibody at 0.5ug/ml. The expected ~43/38 kDa (isoforms beta 1/gamma 1) bands and the additional ~21kDa band are all blocked by the immunizing peptide.



IHC testing of FFPE human spleen tissue with GABPB2 antibody at 3.75ug/ml. Required HIER: steamed antigen retrieval with pH6 citrate buffer; AP-staining.



IHC testing of FFPE human tonsil tissue with GABPB2 antibody at 3.75ug/ml. Required HIER: steamed antigen retrieval with pH6 citrate buffer; AP-staining.

Description

Additional name(s) for this target protein: GA binding protein transcription factor, beta subunit 2

Application Notes

Optimal dilution of the GABPB2 antibody should be determined by the researcher.

1. This GABPB2 antibody will detect the beta and gamma isoforms of the protein.

Immunogen

Amino acids SSENSKATDETG were used as the immunogen for this GABPB2 antibody.

Storage

Aliquot and store the GABPB2 antibody at -20oC.