


## GABRB3 Antibody (R33849)

Catalog No.	Formulation	Size
R33849-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2562
Localization	Cytoplasmic, membrane
Applications	Immunohistochemistry (FFPE) : 2-4ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This GABRB3 antibody is available for research use only.



IHC testing of FFPE human cortex with GABRB3 antibody at 2.5ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

### Description

Additional name(s) for this target protein: Gamma-aminobutyric acid (GABA) A receptor, beta 3

### Application Notes

Optimal dilution of the GABRB3 antibody should be determined by the researcher.

## **Immunogen**

Amino acids EKTAKAKNDRSK were used as the immunogen for this GABRB3 antibody.

## **Storage**

Aliquot and store the GABRB3 antibody at -20oC.