


GABRG2 Antibody (R34865)

Catalog No.	Formulation	Size
R34865-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2566
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This GABRG2 antibody is available for research use only.



Western blot testing of mouse brain lysate with GABRG2 antibody at 2ug/ml. Predicted molecular weight: ~55kDa.

Description

Additional name(s) for this target protein: Gamma-aminobutyric acid (GABA) A receptor, gamma 2

Application Notes

Optimal dilution of the GABRG2 antibody should be determined by the researcher.

Immunogen

Amino acids NNATHLQERDE were used as the immunogen for this GABRG2 antibody.

Storage

Aliquot and store the GABRG2 antibody at -20oC.