


GAD67 Antibody (R35446)

Catalog No.	Formulation	Size
R35446-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Dog, Pig, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2571
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This GAD67 antibody is available for research use only.



Description

Additional name(s) for this target protein: Glutamate decarboxylase 67 kDa isoform, glutamate decarboxylase 1 (brain, 67kDa), GAD1 (67 kDa isoform)

Application Notes

Optimal dilution of the GAD67 antibody should be determined by the researcher.

Immunogen

Amino acids QRTNSLEEKSR were used as the immunogen for this GAD67 antibody.

Storage

Aliquot and store the GAD67 antibody at -20oC.