


Galectin 1 Antibody (R34565)

Catalog No.	Formulation	Size
R34565-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3956
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This Galectin 1 antibody is available for research use only.



Western blot testing of mouse [A] and rat [B] pancreas lysates and NIH3T3 [C] lysate with Galectin 1 antibody at 0.1ug/ml. Predicted/observed molecular weight: ~15kDa.

Description

Additional name(s) for this target protein: Lgals1

Application Notes

Optimal dilution of the Galectin 1 antibody should be determined by the researcher.

Immunogen

Amino acids KPGECLKVRGEVASE were used as the immunogen for this Galectin 1 antibody.

Storage

Aliquot and store the Galectin 1 antibody at -20°C.