


Galectin 3 Antibody (R33294)

Catalog No.	Formulation	Size
R33294-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3958
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Galectin 3 antibody is available for research use only.



Western blot testing of human [A] and mouse [B] colon lysates with Galectin 3 antibody at 0.3ug/ml. Predicted molecular weight ~26kDa.

Description

Additional name(s) for this target protein: Lectin, galactoside-binding, soluble, 3; LGALS3

Application Notes

Optimal dilution of the Galectin 3 antibody should be determined by the researcher.

Immunogen

Amino acids CNTKLDNNWGREERQ were used as the immunogen for this Galectin 3 antibody.

Storage

Aliquot and store the Galectin 3 antibody at -20oC.