


Galectin-3 Antibody (R34764)

Catalog No.	Formulation	Size
R34764-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	16854
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This Galectin-3 antibody is available for research use only.



Western blot testing of mouse heart lysate with Galectin-3 antibody at 0.5ug/ml.
Predicted molecular weight ~26kDa.

Description

Additional name(s) for this target protein: Lectin, galactose binding, soluble 3; Lgals3

Application Notes

Optimal dilution of the Galectin-3 antibody should be determined by the researcher.

Immunogen

Amino acids NANRIVLDFRRGND were used as the immunogen for this Galectin-3 antibody.

Storage

Aliquot and store the Galectin-3 antibody at -20oC.