


GALNS Antibody (R36303)

Catalog No.	Formulation	Size
R36303-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2588
Applications	Western Blot : 0.01-0.03ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This GALNS antibody is available for research use only.



Western blot testing of human bone marrow lysate with GALNS antibody at 0.01ug/ml.
Predicted molecular weight: ~58 kDa.

Description

Additional name(s) for this target protein: Galactosamine (N-acetyl)-6-sulfate sulfatase

Application Notes

Optimal dilution of the GALNS antibody should be determined by the researcher.

Immunogen

Amino acids TTHNLEDHTKLP were used as the immunogen for this GALNS antibody.

Storage

Aliquot and store the GALNS antibody at -20oC.