


GALR1 Antibody (R35768)

Catalog No.	Formulation	Size
R35768-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2587
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This GALR1 antibody is available for research use only.



Western blot testing of rat spinal cord lysate with GALR1 antibody at 0.1ug/ml. Predicted molecular weight ~39kDa.

Description

Additional name(s) for this target protein: Galanin receptor 1, GALNR, GALNR1

Application Notes

Optimal dilution of the GALR1 antibody should be determined by the researcher.

Immunogen

Amino acids RSKPGKPRSTTN were used as the immunogen for this GALR1 antibody.

Storage

Aliquot and store the GALR1 antibody at -20oC.