

GALR2 Antibody (R35823)

Catalog No.	Formulation	Size
R35823-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	8811
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This GALR2 antibody is available for research use only.



Western blot testing of human heart lysate with GALR2 antibody at 0.3ug/ml. Predicted molecular weight ~42kDa.

Description

Additional name(s) for this target protein: Galanin receptor 2, GALNR2

Application Notes

Optimal dilution of the GALR2 antibody should be determined by the researcher.

Immunogen

Amino acids RYPLHSRELRTPRN were used as the immunogen for this GALR2 antibody.

Storage

Aliquot and store the GALR2 antibody at -20oC.