


GALR3 Antibody (R35279)


Catalog No.	Formulation	Size
R35279-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	8484
Applications	Western Blot : 1-3ug/ml Immunohistochemistry (FFPE) : 4-6ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This GALR3 antibody is available for research use only.



Western blot testing of human skin lysate with GALR3 antibody at 2ug/ml. Predicted molecular weight ~39kDa.



IHC staining of FFPE human testis with GALR3 antibody at 3.8ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Galanin Receptor 3

Application Notes

Optimal dilution of the GALR3 antibody should be determined by the researcher.

Immunogen

Amino acids QGPEPREGPVHGGE were used as the immunogen for this GALR3 antibody.

Storage

Aliquot and store the GALR3 antibody at -20oC.