


## GAPDH Antibody [clone BG-7] (RQ5051)

Catalog No.	Formulation	Size
RQ5051	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul


Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)


<b>Availability</b>	1-2 weeks
<b>Species Reactivity</b>	Human, Mouse, Rat
<b>Format</b>	Purified
<b>Host</b>	Rabbit
<b>Clonality</b>	Recombinant Rabbit Monoclonal
<b>Isotype</b>	Rabbit IgG
<b>Clone Name</b>	BG-7
<b>Purity</b>	Affinity purified
<b>UniProt</b>	P04406
<b>Localization</b>	Nuclear, cytoplasmic
<b>Applications</b>	Western Blot : 1:3000-1:10000 Immunohistochemistry (FFPE) : 1:100-1:250 Immunofluorescence : 1:100-1:250
<b>Limitations</b>	This GAPDH antibody is available for research use only.




Western blot testing of 1) rat liver, 2) rat brain, 3) rat RH35, 4) mouse liver, 5) mouse brain and 6) mouse NIH 3T3 cell lysate with GAPDH antibody. Predicted molecular weight ~36 kDa.




Western blot testing of 1) human HeLa, 2) human K562, 3) human MCF7, 4) human A549, 5) human PC-3 and 6) monkey COS7 cell lysate with GAPDH antibody. Predicted molecular weight ~36 kDa.



IHC staining of FFPE human bladder cancer with GAPDH antibody. HIER: boil tissue sections in pH6, 10mM citrate buffer, for 10-20 min and allow to cool before testing.



Immunofluorescent staining of human HeLa cells with GAPDH antibody (green) and DAPI nuclear stain (blue).



Immunofluorescent staining of human HeLa cells with GAPDH antibody (green) and DAPI nuclear stain (blue).

## Description

GAPDH antibody is a widely used reagent for research in cell biology, protein expression studies, and molecular assays. The encoded protein, glyceraldehyde 3 phosphate dehydrogenase (GAPDH), is a glycolytic enzyme that catalyzes the conversion of glyceraldehyde 3 phosphate to 1,3 bisphosphoglycerate. Because of its ubiquitous expression and stable abundance, GAPDH is frequently used as a housekeeping protein and loading control in western blotting and quantitative PCR studies.

GAPDH is expressed in virtually all tissues, reflecting its essential role in energy metabolism. Its steady expression across experimental conditions makes it an ideal reference protein for normalization. Detection with GAPDH antibody ensures accurate interpretation of protein loading and experimental outcomes. For this reason, GAPDH is one of the most cited internal controls in biomedical research.

Beyond glycolysis, GAPDH participates in multiple cellular processes. It is involved in DNA repair, nuclear transport, membrane fusion, and apoptosis. GAPDH translocation to the nucleus has been observed during oxidative stress, where it influences gene expression and cell survival. These non glycolytic functions highlight GAPDH as more than just a metabolic enzyme.

At the structural level, GAPDH is a tetrameric protein with conserved catalytic cysteine residues critical for enzymatic activity. Its high degree of evolutionary conservation across species makes GAPDH antibody applicable in a wide range of models, from bacteria to mammals. This versatility further reinforces its value as a reference protein in experimental biology.

The GAPDH antibody is widely applied in western blotting, immunohistochemistry, immunofluorescence, and ELISA. These applications support detection across tissues, cell lines, and experimental systems. For researchers seeking reliable normalization controls or exploring glycolysis and multifunctional protein roles, the GAPDH antibody is indispensable. NSJ Bioreagents supplies validated antibodies that ensure reproducibility and accuracy in molecular studies.

## **Application Notes**

Optimal dilution of the GAPDH antibody should be determined by the researcher.

## **Immunogen**

A synthetic peptide specific to human GAPDH was used as the immunogen for the GAPDH antibody.

## **Storage**

Store the GAPDH antibody at -20°C.