


GATA1 Antibody (R33371)

Catalog No.	Formulation	Size
R33371-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Dog, Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2623
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This GATA1 antibody is available for research use only.



Western blot testing of human PBMC lysate with GATA1 antibody at 0.3ug/ml. Predicted molecular weight ~43 kDa.

Description

Additional name(s) for this target protein: GATA binding protein 1 (globin transcription factor 1)

Application Notes

Optimal dilution of the GATA1 antibody should be determined by the researcher.

Immunogen

Amino acids DAEAYRHSPVFQ were used as the immunogen for this GATA1 antibody.

Storage

Aliquot and store the GATA1 antibody at -20oC.