


GATA2 Antibody (F53142)

Catalog No.	Formulation	Size
F53142-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F53142-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P23769
Applications	Western Blot : 1:2000
Limitations	This GATA2 antibody is available for research use only.



Western blot testing of GATA2 antibody at 1:2000 dilution + K562 lysate. Predicted molecular weight ~51 kDa.

Description

Transcriptional activator which regulates endothelin-1 gene expression in endothelial cells. Binds to the consensus sequence 5'-AGATAG-3'.

Application Notes

Titration of the GATA2 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

This GATA2 antibody was produced from a rabbit immunized with a KLH conjugated synthetic peptide between 385-419 amino acids from the C-terminal region of human GATA2.

Storage

Aliquot the GATA2 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.