


## GATA5 Antibody (R33856)

Catalog No.	Formulation	Size
R33856-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	140628
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This GATA5 antibody is available for research use only.



Western blot testing of HepG2 cell lysate with GATA5 antibody at 0.3ug/ml. Predicted molecular weight ~41kDa.

## Description

Additional name(s) for this target protein: GATA binding protein 5

## Application Notes

Optimal dilution of the GATA5 antibody should be determined by the researcher.

## Immunogen

Amino acids ELAARPGWAQTATAD were used as the immunogen for this GATA5 antibody.

## **Storage**

Aliquot and store the GATA5 antibody at -20oC.