



GCNT3 Antibody (R35573)

Catalog No.	Formulation	Size
R35573-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9245
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This GCNT3 antibody is available for research use only.





Additional name(s) for this target protein: Glucosaminyl (N-acetyl) transferase 3

Application Notes

Optimal dilution of the GCNT3 antibody should be determined by the researcher.

Immunogen

Amino acids RDTLHLTNKKKD were used as the immunogen for this GCNT3 antibody.

Storage

Aliquot and store the GCNT3 antibody at -20oC.