


GGT1 Antibody (R34849)

Catalog No.	Formulation	Size
R34849-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2678
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This GGT1 antibody is available for research use only.



Western blot testing of human kidney lysate with GGT1 antibody at 0.1ug/ml. Predicted molecular weight ~61kDa.

Description

Additional name(s) for this target protein: Gamma-glutamyltransferase 1

Application Notes

Optimal dilution of the GGT1 antibody should be determined by the researcher.

Immunogen

Amino acids NKRTVIEQQPVLCE were used as the immunogen for this GGT1 antibody.

Storage

Aliquot and store the GGT1 antibody at -20oC.