


## GGT1 Antibody (R35149)

Catalog No.	Formulation	Size
R35149-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2678
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This GGT1 antibody is available for research use only.



Western blot testing of human kidney lysate with GGT1 antibody at 0.3ug/ml. Predicted molecular weight ~61kDa.

### Description

Additional name(s) for this target protein: Gamma-glutamyltransferase 1

### Application Notes

Optimal dilution of the GGT1 antibody should be determined by the researcher.

### Immunogen

Amino acids PSASKEPDNHVYTR were used as the immunogen for this GGT1 antibody.

## Storage

Aliquot and store the GGT1 antibody at -20oC.