

## GGT3P Antibody (Center Region) (F54064)

Catalog No.	Formulation	Size
F54064-0.2ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.2 ml
F54064-0.05ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.05 ml

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Format</b>	Antigen affinity purified
<b>Clonality</b>	Polyclonal (rabbit origin)
<b>Isotype</b>	Rabbit Ig
<b>Purity</b>	Antigen affinity
<b>UniProt</b>	A6NGU5
<b>Applications</b>	Western Blot : 1:2000
<b>Limitations</b>	This GGT3P antibody is available for research use only.



Western blot testing of GGT3P antibody at 1:2000 + human kidney lysate. Predicted molecular weight ~62 kDa.

### Description

Initiates extracellular glutathione (GSH) breakdown; catalyzes the transfer of the glutamyl moiety of glutathione to amino acids and dipeptide acceptors. [UniProt]

### Application Notes

The stated application concentrations are suggested starting points. Titration of the GGT3P antibody may be required due to differences in protocols and secondary/substrate sensitivity.

### Immunogen

A portion of amino acids 176-207 from human GGT3P was used as the immunogen for the GGT3P antibody.

## **Storage**

Aliquot the GGT3P antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.