


Glutathione peroxidase 2 Antibody (R35954)

Catalog No.	Formulation	Size
R35954-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Dog, Pig, Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2877
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Glutathione peroxidase 2 antibody is available for research use only.



Western blot testing of human liver lysate with Glutathione peroxidase 2 antibody at 0.5ug/ml. Predicted molecular weight ~22kDa.

Description

Additional name(s) for this target protein: GPX2, GI-GPx, GSHPX-GI

Application Notes

Optimal dilution of the Glutathione peroxidase 2 antibody should be determined by the researcher.

Immunogen

Amino acids SLDGEKVDFN were used as the immunogen for this Glutathione peroxidase 2 antibody.

Storage

Aliquot and store the Glutathione peroxidase 2 antibody at -20oC.