


Glutathione peroxidase Antibody (R35473)

Catalog No.	Formulation	Size
R35473-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	14775
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This Glutathione peroxidase antibody is available for research use only.



Western blot testing of rat liver lysate with Glutathione peroxidase antibody at 0.3ug/ml.
Predicted molecular weight ~22kDa.

Description

Additional name(s) for this target protein: Gpx, Gpx1

Application Notes

Optimal dilution of the Glutathione peroxidase antibody should be determined by the researcher.

Immunogen

Amino acids QSTVYAFSARP were used as the immunogen for this Glutathione peroxidase antibody.

Storage

Aliquot and store the Glutathione peroxidase antibody at -20oC.