


GNL3 Antibody (R33842)

Catalog No.	Formulation	Size
R33842-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	Q9BVP2
Gene ID	26354
Applications	Western Blot : 2-4ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This GNL3 antibody is available for research use only.



Western blot testing of A) HeLa, B) HepG2, and C) K562 nuclear lysates with GNL3 antibody at 2ug/ml. Predicted molecular weight: ~62/60kDa (isoforms 1/2). The observed ~75kDa band likely represents the acetylated form of the protein.

Description

Additional name(s) for this target protein: Guanine nucleotide binding protein-like 3

Application Notes

Optimal dilution of the GNL3 antibody should be determined by the researcher.

Immunogen

Amino acids EKNNAQSIKRAIKGPH were used as the immunogen for this GNL3 antibody.

Storage

Aliquot and store the GNL3 antibody at -20oC.