


## GNRHR Antibody (R35341)

Catalog No.	Formulation	Size
R35341-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2798
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This GNRHR antibody is available for research use only.



Western blot testing of human peripheral blood mononucleocyte lysate with GNRHR antibody at 0.5ug/ml. Expected molecular weight: 38-60 kDa but may be observed at higher molecular weights due to glycosylation.

## Description

Additional name(s) for this target protein: Gonadotropin-releasing hormone receptor, GRHR, LRHR

## Application Notes

Optimal dilution of the GNRHR antibody should be determined by the researcher.

## Immunogen

Amino acids QKWTQKKEKGK were used as the immunogen for this GNRHR antibody.

## Storage

Aliquot and store the GNRHR antibody at -20oC.