


## GOT2 Antibody (R35988)

| Catalog No.  | Formulation   | Size   |
|--------------|---|--------|
| R35988-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

|                      |   |
|----------------------|---|
| Availability         | 1-3 business days   |
| Species Reactivity   | Human   |
| Predicted Reactivity | Mouse, Rat  |
| Format               | Antigen affinity purified                                 |
| Host                 | Goat  |
| Clonality            | Polyclonal (goat origin)                                  |
| Isotype              | Goat Ig   |
| Purity               | Antigen affinity  |
| Gene ID              | 2806  |
| Applications         | Western Blot : 0.3-1ug/ml<br>ELISA (peptide) LOD : 1:4000 |
| Limitations          | This GOT2 antibody is available for research use only.    |



Western blot testing of human liver lysate with GOT2 antibody at 0.3ug/ml. Predicted molecular weight ~47kDa.

### Description

Additional name(s) for this target protein: Glutamic-oxaloacetic transaminase 2

### Application Notes

Optimal dilution of the GOT2 antibody should be determined by the researcher.

## **Immunogen**

Amino acids SNLKKEGSTHNWQH were used as the immunogen for this GOT2 antibody.

## **Storage**

Aliquot and store the GOT2 antibody at -20oC.