


GPR119 Antibody (R35841)

Catalog No.	Formulation	Size
R35841-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	139760
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This GPR119 antibody is available for research use only.



Western blot testing of human ovary lysate with GPR119 antibody at 0.3ug/ml. Predicted molecular weight ~37kDa.

Description

Additional name(s) for this target protein: G protein-coupled receptor 119

Application Notes

Optimal dilution of the GPR119 antibody should be determined by the researcher.

Immunogen

Amino acids SARNCGPERPRESS were used as the immunogen for this GPR119 antibody.

Storage

Aliquot and store the GPR119 antibody at -20oC.