


GPR120 Antibody (F45258)

Catalog No.	Formulation	Size
F45258-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F45258-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q5NUL3
Applications	Western Blot : 1:1000
Limitations	This GPR120 antibody is available for research use only.



Description

GPR120 is a member of the rhodopsin family of G protein-coupled receptors (GPRs).

Application Notes

Titration of the GPR120 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 231-260 from the human protein was used as the immunogen for this GPR120 antibody.

Storage

Aliquot the GPR120 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.