


GPR83 Antibody (R36164)

Catalog No.	Formulation	Size
R36164-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10888
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This GPR83 antibody is available for research use only.



Western blot testing of human cerebellum lysate with GPR83 antibody at 1ug/ml.
Predicted molecular weight: ~48kDa.

Description

Additional name(s) for this target protein: G protein-coupled receptor 83

Application Notes

Optimal dilution of the GPR83 antibody should be determined by the researcher.

Immunogen

Amino acids TEPHEGRADEQSAE were used as the immunogen for this GPR83 antibody. The amino acid sequence used as immunogen is from the extracellular portion of the protein.

Storage

Aliquot and store the GPR83 antibody at -20oC.