


Gpx1 Antibody (R33689)

Catalog No.	Formulation	Size
R33689-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse, Rat
Predicted Reactivity	Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	14775
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Gpx1 antibody is available for research use only.



Western blot testing of rat liver lysate with Gpx1 antibody at 0.1g/ml. Predicted molecular weight ~22kDa.

Description

Additional name(s) for this target protein: Glutathione peroxidase 1

Application Notes

Optimal dilution of the Gpx1 antibody should be determined by the researcher.

Immunogen

Amino acids DYTEMNDLQKRLGPR were used as the immunogen for this Gpx1 antibody.

Storage

Aliquot and store the Gpx1 antibody at -20oC.