


## Gpx1 Antibody (R33689)

Catalog No.	Formulation	Size
R33689-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse, Rat
Predicted Reactivity	Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	14775
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Gpx1 antibody is available for research use only.



Western blot testing of rat liver lysate with Gpx1 antibody at 0.1g/ml. Predicted molecular weight ~22kDa.

## Description

Additional name(s) for this target protein: Glutathione peroxidase 1

## Application Notes

Optimal dilution of the Gpx1 antibody should be determined by the researcher.

## Immunogen

Amino acids DYTEMNDLQKRLGPR were used as the immunogen for this Gpx1 antibody.

## **Storage**

Aliquot and store the Gpx1 antibody at -20oC.