


## GRAMD3 Antibody (R35490)

Catalog No.	Formulation	Size
R35490-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Dog, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	65983
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This GRAMD3 antibody is available for research use only.



Western blot testing of human temporal cortex lysate with GRAMD3 antibody at 1ug/ml.  
Predicted molecular weight: ~48kDa.

### Description

Additional name(s) for this target protein: GRAM domain containing 3

### Application Notes

Optimal dilution of the GRAMD3 antibody should be determined by the researcher.

## Immunogen

Amino acids QTPESSENSRD were used as the immunogen for this GRAMD3 antibody.

## Storage

Aliquot and store the GRAMD3 antibody at -20oC.