


GRAP Antibody (R34086)

Catalog No.	Formulation	Size
R34086-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10750
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This GRAP antibody is available for research use only.



Western blot testing of human thymus lysate with GRAP antibody at 0.3ug/ml. Predicted molecular weight: ~25kDa.

Description

Additional name(s) for this target protein: GRB2-related adaptor protein

Application Notes

Optimal dilution of the GRAP antibody should be determined by the researcher.

Immunogen

Amino acids GFFPRSYVQPVHL were used as the immunogen for this GRAP antibody.

Storage

Aliquot and store the GRAP antibody at -20oC.