


GRB7 Antibody (R33858)


Catalog No.	Formulation	Size
R33858-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2886
Applications	Western Blot : 2-5ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This GRB7 antibody is available for research use only.



Western blot testing of A431 lysate with GRB7 antibody at 2ug/ml. Predicted molecular weight: ~60kDa.



Western blot of HEK293 lysate overexpressing GRB7 probed with GRB7 antibody (mock transfection in lane 1). Predicted molecular weight: ~60kDa.

Description

Additional name(s) for this target protein: Growth factor receptor-bound protein 7

Application Notes

Optimal dilution of the GRB7 antibody should be determined by the researcher.

1. This GRB7 antibody will detect multiple isoforms from 173-555 amino acids in length.

Immunogen

Amino acids ELDLSPPHLSSSPE were used as the immunogen for this GRB7 antibody.

Storage

Aliquot and store the GRB7 antibody at -20oC.