


GRIN3B Antibody (R33711)

Catalog No.	Formulation	Size
R33711-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	116444
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This GRIN3B antibody is available for research use only.



Western blot testing of human cerebellum lysate with GRIN3B antibody at 0.5ug/ml.
Predicted molecular weight: ~113kDa.

Description

Additional name(s) for this target protein: Glutamate receptor, ionotropic, N-methyl-D-aspartate 3B

Application Notes

Optimal dilution of the GRIN3B antibody should be determined by the researcher.

Immunogen

Amino acids EAPPHSGRPGSQ were used as the immunogen for this GRIN3B antibody.

Storage

Aliquot and store the GRIN3B antibody at -20oC.