


## GSN Antibody (R35137)

Catalog No.	Formulation	Size
R35137-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2934
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This GSN antibody is available for research use only.



Western blot testing of human spleen lysate with GSN antibody at 0.1ug/ml.  
Predicted/observed molecular weight ~86kDa.

## Description

Additional name(s) for this target protein: Gelsolin

## Application Notes

Optimal dilution of the GSN antibody should be determined by the researcher.

## **Immunogen**

Amino acids PRLKDKKMDAHP were used as the immunogen for this GSN antibody.

## **Storage**

Aliquot and store the GSN antibody at -20oC.