


GSN Antibody (R35137)

Catalog No.	Formulation	Size
R35137-100UG	0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2934
Applications	Western Blot: 0.1-0.3ug/ml ELISA (peptide) LOD: 1:128000
Limitations	This GSN antibody is available for research use only.



Description

Additional name(s) for this target protein: Gelsolin

Application Notes

Optimal dilution of the GSN antibody should be determined by the researcher.

Immunogen

Amino acids PRLKDKKMDAHP were used as the immunogen for this GSN antibody. **Storage** Aliquot and store the GSN antibody at -20oC.