

GSTA3 Antibody (R36490)

Catalog No.	Formulation	Size
R36490-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	Q16772
Applications	Western Blot : 0.2-0.5ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This GSTA3 antibody is available for research use only.



Western blot testing of human placental lysate with GSTA3 antibody at 0.3ug/ml.
Predicted molecular weight ~25 kDa.

Description

Alternate name(s) for this target: GST class-alpha member 3; Glutathione S-transferase A3

Application Notes

Optimal dilution of the GSTA3 antibody should be determined by the researcher.

Immunogen

Amino acids RKPPADAKALE were used as the immunogen for this GSTA3 antibody.

Storage

Aliquot and store the GSTA3 antibody at -20oC.