


GSTA4 Antibody (R36449)

Catalog No.	Formulation	Size
R36449-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	O15217
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) : 1:64000
Limitations	This GSTA4 antibody is available for research use only.



Western blot testing of human liver lysate with GSTA4 antibody at 0.3ug/ml. Predicted molecular weight ~26 kDa.

Description

Additional name(s) for this target protein: Glutathione S-transferase alpha 4

Application Notes

The stated application concentrations are suggested starting points. Titration of the GSTA4 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

Amino acids QEYTVKLSNIPTIK were used as the immunogen for this GSTA4 antibody.

Storage

Aliquot and store the GSTA4 antibody at -20oC.