


Gsta5 Antibody (R33881)

Catalog No.	Formulation	Size
R33881-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse, Rat
Predicted Reactivity	Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	14859
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Gsta5 antibody is available for research use only.



Western blot testing of rat liver lysate with Gsta5 antibody at 0.3ug/ml. Predicted molecular weight: ~25kDa.

Description

Additional name(s) for this target protein: Glutathione S-transferase alpha-5

Application Notes

Optimal dilution of the Gsta5 antibody should be determined by the researcher.

1. This Gsta5 antibody will also detect the homologous mouse Gsta3.

Immunogen

Amino acids EMDPGIVDNFLLK were used as the immunogen for this Gsta5 antibody.

Storage

Aliquot and store the Gsta5 antibody at -20oC.