


GSTM1 Antibody (R33361)

Catalog No.	Formulation	Size
R33361-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2944; 2946
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:1000
Limitations	This GSTM1 antibody is available for research use only.



Western blot testing of human lung lysate with GSTM1 antibody at 0.3ug/ml. Predicted molecular weight: ~26kDa.

Description

Additional name(s) for this target protein: Glutathione S-transferase, Mu; GSTM

Application Notes

Optimal dilution of the GSTM1 antibody should be determined by the researcher.

1. This GSTM1 antibody will detect both GSTM1 and GSTM2.

Immunogen

Amino acids CVFSKMAVWGNK were used as the immunogen for this GSTM1 antibody.

Storage

Aliquot and store the GSTM1 antibody at -20oC.