


GSTO1 Antibody (R34395)

Catalog No.	Formulation	Size
R34395-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9446
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This GSTO1 antibody is available for research use only.



Western blot testing of human liver [A] and peripheral blood lymphocyte [B] lysates with GSTO1 antibody at 0.5ug/ml. Predicted molecular weight: ~28kDa.

Description

Additional name(s) for this target protein: Glutathione S-transferase omega 1

Application Notes

Optimal dilution of the GSTO1 antibody should be determined by the researcher.

Immunogen

Amino acids KEDPTVSALLTSEKD were used as the immunogen for this GSTO1 antibody.

Storage

Aliquot and store the GSTO1 antibody at -20oC.