


H pylori NAP Antibody / Helicobacter pylori Neutrophil-activating protein [clone RM414] (R20430)

Catalog No.	Formulation	Size
R20430-0.1ML	1 mg/ml in PBS with 50% glycerol, 1% BSA and 0.09% sodium azide	100 ul

Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Helicobacter Pylori
Format	Purified
Host	Rabbit
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	RM414
Purity	Protein A purified from animal origin-free supernatant
Applications	ELISA : 0.01ug/ml-0.2ug/ml
Limitations	This recombinant H pylori NAP antibody is available for research use only.



An ELISA of Helicobacter Pylori proteins using recombinant H pylori NAP antibody. The plate was coated with 1 ug/ml of CagA, OMP, Urease, or HP-NAP of H. Pylori. A serial dilution of the RM414 mAb was used as the primary antibody. An alkaline phosphatase conjugated anti-rabbit IgG was used as the secondary antibody.

Description

This antibody reacts to the neutrophil-activating protein of Helicobacter pylori.

Application Notes

The stated application concentrations are suggested starting points. Titration of the recombinant H pylori NAP antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

Proteins purified from *Helicobacter pylori* was used as the immunogen for the recombinant H pylori NAP antibody.

Storage

Store the recombinant H pylori NAP antibody at -20oC.