


H1L Antibody [clone 20CT26.44.10.3] (F40184)

Catalog No.	Formulation	Size
F40184-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F40184-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Vaccinia Virus
Format	Purified
Host	Mouse
Clonality	Monoclonal (mouse origin)
Isotype	Mouse IgG2a
Clone Name	20CT26.44.10.3
Purity	Purified
UniProt	P20495
Applications	Western Blot : 1:100-1:500
Limitations	This H1L antibody is available for research use only.



Description

This Protein shows activity toward tyrosine-protein phosphate as well as with serine-protein phosphate.

Application Notes

Titration of the H1L antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

Purified recombinant H1L fusion protein was used to produce this monoclonal antibody.

Storage

Aliquot the H1L antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.