


H1L Antibody [clone 20CT26.44.10.3] (F40184)

Catalog No.	Formulation	Size
F40184-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F40184-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

Bulk quote request

Availability	1-3 business days
Species Reactivity	Vaccinia Virus
Format	Purified
Clonality	Monoclonal (mouse origin)
Isotype	Mouse IgG2a
Clone Name	20CT26.44.10.3
Purity	Purified
UniProt	P20495
Applications	Western Blot : 1:100-1:500
Limitations	This H1L antibody is available for research use only.



Description

This Protein shows activity toward tyrosine-protein phosphate as well as with serine-protein phosphate.

Application Notes

Titration of the H1L antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

Purified recombinant H1L fusion protein was used to produced this monoclonal antibody. **Storage** Aliquot the H1L antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles. Ordering: Phone: 858.663.9055 | Fax: 1.267.821.0800 | Email: info@nsjbio.com Copyright $\ensuremath{\texttt{@}}$ NSJ Bioreagents. All rights reserved.